


产品名称: GMP Recombinant Anti-Human CD3 Monoclonal Antibody
产品货号: AMRe90003



概述 (Summary)

英文全称	GMP Recombinant Anti-Human CD3 Monoclonal Antibody
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤ 10 EU/mg
抗体亚型 (Isotype)	IgG2a/kappa is expressed.
表达宿主 (Host)	Mammalian cell
预测分子量 (Predicted MW)	150 kDa
产品形态 (Form)	Liquid
储存缓冲液 (Buffer)	0.22 μ m-filtered solution containing 20 mM sodium phosphate, 150 mM sodium chloride, pH 7.4
运输方式 (Shipping)	The product is shipped with blue ice.
稳定性&储存 (Stability&Storage)	6 months at 2°C to 8°C under sterile conditions 24 months at -20°C to -80°C under sterile conditions Use a manual defrost freezer and avoid repeated freeze-thaw cycles.

电泳图 (SDS-PAGE image)



注意事项 (Note)

For research use only .