

产品名称: Ra-blot Anti-Rabbit IgG for IP (HRP)
产品货号: AMS0060



产品概述 (Summary)

产品名称 (Production Name)	Ra-blot Anti-Rabbit IgG for IP (HRP)
应用 (Application)	WB,IP,ELISA
种属反应性 (Reactivity)	Rabbit
偶联物 (Conjugation)	HRP
同种型 (Isotype)	IgG2b kappa
克隆 (Clonality)	Monoclonal Antibody
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, 50% glycerol, pH 7.4.

免疫原信息 (Immunogen)

The antibody was developed in goat using the nature rabbit IgG as the immunogen and then develop monoclonal antibody.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2500-1:10000 , ELISA 1:1000-1:3000
------------------------	---

背景信息 (Background)

Ra-blot Anti-Rabbit IgG secondary antibody is western blotting reagents that enable the trouble-free detection of western blotted target protein bands from upstream IP/CoIP, without interference from denatured IgG (IP primary antibody). This allows to detect the (co-)immunoprecipitated protein without interfering by the IgG heavy (50 kDa) and light chains (25 kDa).