

产品名称: Fluor405-conjugated Polyclonal Goat Anti-Mouse IgG(H+L)

产品货号: APS0101

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

Production Name	Fluor405-conjugated Polyclonal Goat Anti-Mouse IgG(H+L)
Clonality	Polyclonal
Host	Goat
Reactivity	Mouse
Application	IF, FC
Form	Liquid
Concentration	0.5mg/mL
Immunogen	IgG
Specificity	IgG(H+L)
Buffer	50% Glycerol,0.01M PBS,0.5%BSA,0.1%PC300,pH7.4
Conjugation	Fluor405
Ex-max	400nm
Em-max	424nm
Purification	Affinity purification
Purity	Calculated as antibodies, no less than 90%.
Storage	Please store at -20℃ away from light.

使用建议

IF:1:20-100,FC:1~10ug/ml,This data is a recommended range, please refer to the actual value according to the specific situation.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.