

产品名称: PRKG1 (cGKI) (7A16) Mouse Monoclonal antibody

产品货号: AMM25628

产品概述 (Summary)

产品名称 (Production Name)	PRKG1 (cGKI) (7A16) Mouse Monoclonal antibody
描述 (Description)	Mouse Monoclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	IF-P,IF-F,IF-ICC,IHC-P,IHC-F,WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG1
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% New type preservative N.
纯化方式 (Purification)	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)

免疫原信息 (Immunogen)

基因名 (Gene Name)	cGKI
别名 (Alternative Names)	cGKI

产品名称: PRKG1 (cGKI) (7A16) Mouse Monoclonal antibody

产品货号: AMM25628

蛋白 ID (SwissProt ID)

NP_006249

产品应用 (Application)


稀释比 (Dilution Ratio)

WB 1:3000, IHC-P/IHC-F 1:50, IF-P/IF-F/IF-ICC 1:100


蛋白分子量 (Molecular Weight)

Observed MW:77.6kDa

图片信息 (Image)




Immunohistochemical staining of paraffin-embedded Kidney tissue within the normal limits using anti-PRKG1 mouse monoclonal antibody.




Immunohistochemical staining of paraffin-embedded pancreas tissue within the normal limits using anti-PRKG1 mouse monoclonal antibody.

产品名称: PRKG1 (cGKI) (7A16) Mouse Monoclonal antibody

产品货号: AMM25628



Anti-PRKG1 mouse monoclonal antibody immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PRKG1



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PRKG1 (Right lane) cDNA for 48 hrs and lysed.

注意事项 (Note)

For research use only .