

产品名称: STING (Phospho Ser366) Rabbit Polyclonal Antibody
产品货号: AP Rab59196



产品概述 (Summary)

产品名称 (Production Name)	STING (Phospho Ser366) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit Polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat,Mouse,

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Phospho
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
纯化方式 (Purification)	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMEM173
别名 (Alternative Names)	ERIS MITA STING
基因 ID (Gene ID)	Human:340061

产品名称: STING (Phospho Ser366) Rabbit Polyclonal Antibody
产品货号: AP Rab59196

蛋白 ID (SwissProt ID)

Human:Q86WV6,Mouse:Q3TBT3.Synthesized phospho peptide around human
STING (Ser366)

产品应用 (Application)

稀释比 (Dilution Ratio)

WB 1:1000-2000

蛋白分子量 (Molecular Weight)

Observed MW:38kDa


研究背景 (Background)

This gene encodes a five transmembrane protein that functions as a major regulator of the innate immune response to viral and bacterial infections. The encoded protein is a pattern recognition receptor that detects cytosolic nucleic acids and transmits signals that activate type I interferon responses. The encoded protein has also been shown to play a role in apoptotic signaling by associating with type II major histocompatibility complex. Mutations in this gene are the cause of infantile-onset STING-associated vasculopathy. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2014],

研究领域 (Research Area)

-

图片 (Image Date)



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .